

FOR IMMEDIATE RELEASE

Nov. 11, 2009

Totoku Launches the CCL240, its New, Color, 24.1", 2.3MP, Widescreen LCD for Clinical Review

November 09 2009 – Tokyo, Japan - Totoku is pleased to announce its new CCL240, a 2.3MP, wide-screen, DICOM LCD, designed as a clinical review system for radiology workstations and beyond.

"We believe the CCL240 can be placed in virtually any department where secondary review may occur," says Mark Yoshida, Product Manager of Totoku North America. "Departments such as cardiology, oncology, nuclear medicine, and pathology will appreciate the wide viewing angles of the display, which allow physicians to view multiple images on a single display, simultaneously."



The CCL240's ultra-wide 24.1" Super-IPS panel helps prevent color shift at the furthest vantage points, thereby improving viewing angles for impeccable image reproduction and a superior ergonomic design. Additionally, the ultra-wide panel enables physicians to view multiple images from multiple modalities, simultaneously.

The CCL240 is compatible with any PACS, cardiology, EMR, or pathology network and is DICOM calibrated for confident image review.

Additional Specifications:

- 1920 x 1200 resolution
- Full DICOM calibration
- 10-bit built-in LUT; luminance stabilization system
- S-IPS panel with ultra-wide viewing angles
- 400cd/m² brightness
- 1000:1 contrast ratio
- Optional A/R coating for sharp, film-like black hues

The CCL240 will begin shipping December 2009, will be on exhibit at RSNA 2009 in Chicago at the following booths:

- Totoku / US Electronics – booth #4065
- Quest International, Inc. – booth #7303

About TOTOKU

TOTOKU, founded in 1940 and headquartered in Tokyo, Japan, is a leading manufacturer of high-end displays for medical imaging and industrial markets worldwide. Since its launch on the display business in 1972, TOTOKU has supplied its specialty displays into multiple vertical markets and has earned high marks for its product reliability and post installation support, in addition to specialized technologies required for the respective fields. TOTOKU's products and services are available via its US, EU, and Asia offices and distributors throughout the world.

USA	TOTOKU NORTH AMERICA INC., in USA	TEL: +1-469-948-4839
EUROPE	TOTOKU Business Unit. in Germany	TEL: +49-2156-496880
ASIA	TOTOKU ELECTRIC CO., LTD. in Japan	TEL: +81-3-5860-2132

**###
 END**