

FOR IMMEDIATE RELEASE

June 10, 2009

## Totoku Launches CCL208, a Feature-Rich, Color LCD for CT and MR Applications

TOTOKU offers the CCL208, a new 20,1", color display, with a resolution of two megapixels, a brightness of 300cd/m<sup>2</sup> and a contrast ratio of 800: 1 it's suitable for all CT and MR applications as well as dental x-ray images.

Noteworthy features include: integrated 10-Bit LUT (look up table) for precise color and grayscale reproduction and full hardware calibration, the most precise calibration method in the market today. "The CCL208 includes Medivisor™ Color calibration software which offers a full DICOM conformance check and calibration", explains Takahiro Takayama, General manager of Medical Displays Sales and Marketing division, Totoku Electric Company. "The feature-rich CCL208 achieves a perfect balance of cost and performance."

Equipped with the Lambda-Sentinel-I circuit, all brightness differences due to aging or temperature changes are compensated. Employing this technology, after thirty seconds, the display reaches full diagnostic brightness level. In addition, high diagnostic precision is secured over the complete lifetime of the display.

In addition to image quality, the CCL208 demonstrates superior price performance providing the best value for a preset DICOM display, with full hardware calibration in the market today.

### About TOTOKU

TOTOKU, founded in 1940 and headquartered in Tokyo, Japan, is a leading manufacturer of high-end displays for medical imaging and industrial markets worldwide. Since its launch on the display business in 1972, TOTOKU has supplied its specialty displays into multiple vertical markets and has earned high marks for its product reliability and post installation support, in addition to specialized technologies required for the respective fields. TOTOKU's products and services are available via its US, EU, and Asia offices and distributors throughout the world.

USA	TOTOKU NORTH AMERICA INC., in USA	TEL: +1-972-550-1888
EUROPE	TOTOKU Europe GmbH. in Germany	TEL: +49-2131-36668-0
ASIA	TOTOKU ELECTRIC CO., LTD. in Japan	TEL: +81-3-5273-2053

###  
END